

Product	Size (cm)	BUILD HEIGHT ¹		WEIGHT ³	
		Inches	Millimeters	Ounces	Grams
Sierra (FS1)	22-25	5 1/2	140	11	320
	26-28	6	154	14	390
	29-31	6 3/4	172	16	440
Highlander (FS3)	22-25	3 7/8	98	10	275
	26-28	4 1/4	108	12	340
	29-31	4 3/4	121	16	440
Pacifica (FS2)	22-25	2 3/8	59	10	290
	26-28	2 3/8	61	13	360
	29-31	2 1/2	63	14	410
Pacifica LP (FS4)	22-25	1 7/8	47	10	270
	26-28	2	50	11	320
	29-31	2 3/8	60	16	450
Renegade (RS1)	22-25	6 1/8	154	15	430
	26-28	6 3/4	172	18	515
	29-31	7 1/2	189	22	635
Renegade LP (RS6)	22-25	4 3/4	120	15	430
	26-28	5 1/8	131	18	510
	29-31	5 5/8	143	22	620
Renegade MX (RS1-M1)	22-25	6 5/8	170	29	825
	26-28	7 3/8	188	33	945
	29-31	8 1/8	206	35	990
Renegade LP-MX (RS6-M1)	22-25	5 5/8	143	27	775
	26-28	6	153	30	850
	29-31	6 1/2	166	34	965
Renegade SX (RS1-S1)	22-25	7 1/8	182	27	770
	26-28	7 7/8	200	31	870
	29-31	8 5/8	218	34	965
Renegade LP-SX (RS6-S1)	22-25	6	151	25	700
	26-28	6 1/4	159	27	770
	29-31	6 3/4	172	32	920
Runway/ HX (RS2 & RS2-K1)	22-25	3 5/8	92	15	425
	26-28	3 7/8	98	17	485
Silhouette (RS4 & RS8) w/female adapter	22-25	7	178	12	330
	26-28	7 1/2	191	14	395
	29-31	7 7/8	201	17	475
Silhouette LP (RS5 & RS9) w/male adapter	22-25	6 5/8	168	8	235
	26-28	7 1/8	181	11	300
	29-31	7 1/2	191	13	375
Silhouette VS (R14 & R18) w/female adapter	22-25	7	178	12	330
	26-28	7 1/2	191	14	395
	29-31	7 7/8	201	17	475
Silhouette LP-VS (R15 & R19) w/male adapter	22-25	6 5/8	168	8	235
	26-28	7 1/8	181	11	300
	29-31	7 1/2	191	13	375
Promenade (VS2)	22-25	4 5/8	116	10	295
	26-28	4 7/8	124	12	345
	29-31	5 3/8	135	14	400
LP Symes (LP2) ²	22-25	1 3/4	43	8	225
	26-28	1 7/8	48	10	295
	29-31	2 1/8	53	13	360
Senator (VS1)	22-25	4 3/4	120	11	325
	26-28	5 3/8	136	14	405
	29-31	5 7/8	149	17	470
Defender (FS3)	16-18cm	2 3/8	61	3	90
	19-21cm	3	76	4	105
Freestyle	22-25	4	102	17	485
	26-28	4 1/4	107	19	540
	29-31	4 5/8	117	22	635
Slalom	custom	5 1/2	140	37	1055
Nitro	Tall (FS6)	7 1/4	184	10	275
	Short (FS9)	9 1/2	241	11	300
Foot Shell w/sock (average weights)	22-25		0	7	195
	26-28		0	9	250
	29-31		0	11	310

¹ Measurements are representative of fully assembled units including foot shell and protective sock(s) if applicable. Build height is measured from the bottom of the heel to the middle of the pyramid skirt for male adapters, or to the top of female adapters. Build height will vary based on category selected.

² The LP Symes foot is measured to the top of the posterior surface of the lamination cup.

³ Actual weight varied based on category selected. The reported measurements are average weights.