

2008 Component List

| <u>Part Number</u> | <u>Description</u> | <u>Weight Rating (Lbs.)</u> | <u>Weight Rating (Kg)</u> |
|--------------------|---|-----------------------------|---------------------------|
| ACC-00-10112-00 | Spectra Sock 12" x 2.75" | - | - |
| ACC-00-10114-00 | Spectra Sock 14" x 2.75" | - | - |
| ACC-00-10118-00 | Spectra Sock 18" x 2.75" | - | - |
| KIT-00-10200-00 | Dacron Strap Replacements, 1 st Generation Freedom Series Feet | - | - |
| ACC-00-10300-00 | Foot Shell Removal Tool | - | - |
| ACC-00-10400-00 | Runway Shim Adjustment Wrench | - | - |
| ACC-00-10500-RL | Runway Ankle Mechanism, Replacement, Left | 255 lbs. | 116 kg |
| ACC-00-10500-RR | Runway Ankle Mechanism, Replacement, Right | 255 lbs. | 116 kg |
| RS1-SP | Renegade Sole Plate, Replacement | 365 lbs. | 166 kg |
| RS4-SP | Silhouette Sole Plate, Replacement | 365 lbs. | 166 kg |
| ACC-00-12300-00 | Tube Clamp, 34 mm, Female Pyramid, for Silhouette | 365 lbs. | 165 kg |
| ACC-EX-12300-00 | Extended Tube Clamp, 34mm, Female Pyramid for Silhouette | 365 lbs. | 165 kg |
| KIT-00-123S0-00 | Spacer Ring Set for use with Versa Clamp, Plastic | 365 lbs. | 165 kg |
| KIT-00-123V0-00 | Versa Clamp Adapter Kit, for use with Silhouette | 365 lbs. | 165 kg |
| ACC-00-12400-00 | Insert to 34mm Tube Clamp, Silhouette, Custom Fabrication | 365 lbs. | 165 kg |
| ACC-00-12500-00 | Tube Clamp, 30 mm, Female Pyramid | 300 lbs. | 136 kg |
| ACC-00-12600-00 | Tube Clamp, 30 mm, Male Pyramid | 300 lbs. | 136 kg |
| ACC-00-12700-00 | Tube w/Female Pyramid, 30 mm x 100 mm Titanium | 300 lbs. | 136 kg |
| ACC-00-12800-00 | Tube w/ Titanium Female Pyramid, 30 mm x 400 mm Aluminum | 300 lbs. | 136 kg |
| ACC-00-13000-00 | Pyramid, 4 Hole, Male | 300 lbs. | 136 kg |
| ACC-00-13200-00 | Pyramid, 4 Hole, Female | 300 lbs. | 136 kg |
| ACC-00-13500-00 | Male Pyramid, Rotating, 4 Hole Base | 300 lbs. | 136 kg |
| ACC-00-13600-00 | Female Pyramid, Rotating, 4 Hole Base | 300 lbs. | 136 kg |
| ACC-00-13700-00 | Base, 4 Hole, Rotating Pyramid | 300 lbs. | 136 kg |
| ACC-00-13800-00 | Base, 3 Prong, Rotating Pyramid, Stainless Steel | 220 lbs. | 99 kg |
| ACC-00-13900-00 | Pyramid, Rotating, Male | 300 lbs. | 136 kg |
| ACC-00-14000-00 | Pyramid, Rotating, Female | 300 lbs. | 136 kg |
| ACC-00-14500-00 | Tube, 30 mm graphite, 15", 3 mm wall, Aluminum Inserts | 220 lbs. | 99 kg |
| ACC-00-14600-00 | Tube, 30 mm graphite, 9", 3 mm wall, Aluminum Inserts | 220 lbs. | 99 kg |
| ACC-00-14815-00 | Tube, 30 mm graphite, 15", 4 mm wall, Aluminum Inserts | 300 lbs. | 136 kg |
| ACC-00-14809-00 | Tube, 30 mm graphite, 9", 4 mm wall, Aluminum Inserts | 300 lbs. | 136 kg |
| ACC-00-15000-00 | Pyramid, 4 Prong, Male | 300 lbs. | 136 kg |
| ACC-00-16000-00 | Lamination Ring, 4 Hole | 300 lbs. | 136 kg |
| ACC-00-16100-00 | LP2 Lamination Adaptor, Carbon Fiber (w/ two 10mm bolts) | 365 lbs. | 165 kg |
| ACC-00-17032-00 | Adapter, Double Female Pyramid, 32 mm | 300 lbs. | 136 kg |
| ACC-00-17045-00 | Adapter, Double Female Pyramid, 45 mm | 300 lbs. | 136 kg |
| ACC-00-17060-00 | Adapter, Double Female Pyramid, 60 mm | 300 lbs. | 136 kg |
| ACC-00-17075-00 | Adapter, Double Female Pyramid, 75 mm | 300 lbs. | 136 kg |
| KIT-00-18000-00 | Assembly, Male Pyramid, Cat. 1-6, (Pyramid, 10mm Bolt and Washer) | 220 lbs. | 99 kg |
| KIT-00-18100-00 | Assembly, Male Pyramid, Cat. 7-9, (Pyramid, 12mm Bolt and Washer) | 365 lbs. | 165 kg |
| KIT-00-18200-00 | Assembly, Female Pyramid, Cat.1-6, (Pyramid Receiver, 10mm Bolt, Nut, Washer and 4 Set Screws) | 220 lbs. | 99 kg |
| KIT-00-18300-00 | Assembly, Female Pyramid, Cat.7-9, (Pyramid Receiver, 12mm Bolt, Nut, Washer and 4 Set Screws) | 365 lbs. | 165 kg |
| ACC-00-19500-00 | Knee Pyramid Adapter, Male 30mm , Aluminum | 275 lbs. | 124 kg |
| ACC-00-19600-00 | Male Pyramid Adapter, for use with Silhouette LP | 365 lbs. | 165 kg |
| ACC-00-20100-00 | Alignment Jig Pyramid, Assembly with 10mm bolt (for use with category 6-9 foot products), Stainless Steel | N/A | N/A |
| KIT-00-20200-00 | Heavy Duty Kit, Titanium (incl. two 202-CLA, one 202-EAR, one 202-PYL) | 500 lbs. | 226 kg |
| AFP-00-FS1 | Pre-cut Shin Finishing Foam, for Sierra | - | - |
| AFP-00-RS1 | Pre-cut Shin Finishing Foam, for Renegade (RS 1000 only) | - | - |

* Adapters are constructed of titanium unless noted otherwise.