

# Foot Product Build Height & Weight

Product	Size (cm)	Build Height <sup>1</sup>		Weight <sup>2</sup>	
		Inches	Millimeters	Ounces	Grams
Sierra (FS1)	22-25	5 1/2	140	18	515
	26-28	6	154	23	640
	29-31	6 3/4	172	26	750
Highlander (FS3)	22-25	3 7/8	98	17	470
	26-28	4 1/4	108	21	590
	29-31	4 3/4	121	26	750
Pacifica (FS2)	22-25	2 3/8	59	17	485
	26-28	2 3/8	61	22	610
	29-31	2 1/2	63	25	720
Pacifica LP (FS4)	22-25	1 7/8	47	16	465
	26-28	2	50	20	570
	29-31	2 3/8	60	27	760
Renegade (RS1) & Renegade A-T (R11)	22-25	6 1/8	154	22	625
	26-28	6 3/4	172	27	765
	29-31	7 1/2	189	33	945
Renegade LP (RS6) & Rengade LP-A-T (R16)	22-25	4 3/4	120	22	625
	26-28	5 1/8	131	27	760
	29-31	5 5/8	143	33	930
Renegade MX (RS1-M1)	22-25	6 5/8	170	36	1020
	26-28	7 3/8	188	42	1195
	29-31	8 1/8	206	46	1300
Renegade LP-MX (RS6-M1)	22-25	5 5/8	143	34	970
	26-28	6	153	39	1100
	29-31	6 1/2	166	45	1275
Thrive VS <sup>2</sup> (FS5)	22-25	5 3/4	145	22	630
	26-28	6 1/5	157	27	760
	29-31	7	175	31	870
Runway/ HX (RS2 & RS2-K1)	22-25	3 5/8	92	22	620
	26-28	3 7/8	98	26	735
Silhouette (RS4 & RS8) w/female adapter	22-25	7	178	19	525
	26-28	7 1/2	191	23	645
	29-31	7 7/8	201	28	785
Silhouette LP (RS5 & RS9) w/male adapter	22-25	6 5/8	168	15	430
	26-28	7 1/8	181	19	550
	29-31	7 1/2	191	24	685
Silhouette VS (R14 & R18) w/female adapter	22-25	7	178	19	525
	26-28	7 1/2	191	23	645
	29-31	7 7/8	201	28	785
Silhouette LP-VS (R15 & R19) w/male adapter	22-25	6 5/8	168	15	430
	26-28	7 1/8	181	19	550
	29-31	7 1/2	191	24	685
Promenade (VS2)	22-25	4 5/8	116	17	490
	26-28	4 7/8	124	21	595
	29-31	5 3/8	135	25	710
LP Symes (LP2) <sup>3</sup>	22-25	1 3/4	43	15	420
	26-28	1 7/8	48	19	545
	29-31	2 1/8	53	24	670
Senator (VS1)	22-25	4 3/4	120	18	520
	26-28	5 3/8	136	23	655
	29-30	5 7/8	149	28	780
WalkTek (VS3)	21-25	3 3/8	85	17	490
	26-30	3 3/4	97	22	630
Defender (FS3)	16-18	2 3/8	61	6	175
	19-21	3	76	8	225
Freestyle Swim	22-25	4	102	24	680
	26-28	4 1/4	107	28	790
	29-31	4 5/8	117	33	945
Slalom Ski	custom	5 1/2	140	37	1055
Nitro Running	Tall (FS6)	7 1/4	184	10	275
	Short (FS9)	9 1/2	241	11	300
Foot Shell w/sock (average weights)	22-25			7	195
	26-28			9	250
	29-31			11	310

<sup>1</sup> Measurements are representative of fully assembled units including foot shell and protective sock(s) if applicable. Build height is measured from the bottom of the heel to the middle of the pyramid skirt for male adapters, or to the top of female adapters. Build height will vary based on category selected.

<sup>2</sup> Actual weight varies based on category selected. The reported measurements are average weights.

<sup>3</sup> The LP Symes foot is measured to the top of the posterior surface of the lamination cup.