

PACIFICA® PROSTHETIC FOOT (FS2)

The perfect balance of flexibility and durability

- Designed for amputees with long residual limbs or those needing more build space
- Split toe design provides stability enhancing inversion and eversion
- Lightweight and durable
- Foot lies within the foot shell, so cosmetic finishing is simplified
- Choose from nine stiffness categories and Low, Moderate, and High Impact Levels
- Rated up to 365 lbs



PACIFICA® (FS2)

PACIFICA® PROSTHETIC FOOT (FS2)

Size*: 24-30 cm

*Special Order Size 31 cm (7-10 day lead time)


User Weight Rating: 365 lbs

Stiffness Category: 1-9

Heel Height: 3/8"

Connector: Integral male pyramid

Warranty: 36 months



Clearance & Weight:

Size	Clearance	Avg. Weights	
		Foot	Shell
24-25 cm	2.1 in	259g	127g
26-28 cm	2.2 in	310g	159g
29-31 cm	2.4 in	382g	219g

Note: A Spectra sock and pyramid connector are included with purchase.

REQUIRED ACCESSORIES

Cosmetic Foot Shell:

Mid-profile (Default) FTC-2M

High-profile (Optional) FTC-2H

Choice of caucasian or ebony skin tone

OPTIONAL ACCESSORIES

- **Foot Shell Removal Tool-** Part #ACC-00-103
- **Additional Spectra Sock-** Part #ACC-00-101xx
(Choose 12", 14", or 18" length)
- **Additional Stiffening Bumpers-**
 - KIT-00-111 40 Durometer, Sizes 24-25cm
 - KIT-00-113 40 Durometer, Sizes 26-28cm
 - KIT-00-115 40 Durometer, Sizes 29-31cm