

## PACIFICA® LP PROSTHETIC FOOT (FS4)

*The proven performance of the Pacifica (FS2) in a lower profile design*

- Requires less than 2" of clearance for most commonly used sizes
- Split toe design provides stability enhancing inversion and eversion
- Lightweight and durable
- Foot lies within the foot shell, so cosmetic finishing is simplified
- Choose from nine stiffness categories and Low, Moderate, and High Impact Levels
- Rated up to 365 lbs



PACIFICA® LP (FS4)

# PACIFICA® LP PROSTHETIC FOOT (FS4)

**Size\*: 22-30 cm**

\*Special Order Size 31 cm (7-10 day lead time)


**User Weight Rating:** 365 lbs

**Stiffness Category:** 1-9

**Connector:** Integral male pyramid

**Heel Height:** 3/8"

**Warranty:** 36 months



## Clearance & Weight:

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	1.8 in	305g	127g
26-28 cm	2.0 in	370g	159g
29-31 cm	2.2 in	435g	219g

Note: A Spectra sock and pyramid connector are included with purchase.

## REQUIRED ACCESSORIES

### Cosmetic Foot Shell:

Mid-profile (Default) FTC-2M

High-profile (Optional) FTC-2H

Choice of caucasian or ebony skin tone

## OPTIONAL ACCESSORIES

- **Foot Shell Removal Tool- Part #ACC-00-103**

- **Additional Spectra Sock- Part #ACC-00-101xx**  
(Choose 12", 14", or 18" length)

- **Additional Stiffening Bumpers-**

KIT-00-111 40 Durometer, Sizes 22-25cm

KIT-00-113 40 Durometer, Sizes 26-28cm

KIT-00-115 40 Durometer, Sizes 29-31cm