

Microprocessor-Controlled Knee

Product	Key Features	Avg. Clearance	Avg. Weight	Weight Rating	Alignment	Warranty
Plié® 2.0 MPC Knee	Advanced stumble recovery, water resistant, user replaceable battery	9.25 in / 235 mm	1,235 g	275 lb / 125 kg high impact: 220 lb / 100 kg	Stable, 0-5 mm anterior	36 mo (12 & 24 mo service)

Carbon Fiber Prosthetic Feet

Product	Key Features	Avg. Clearance*	Avg. Weight	Weight Rating	Sizes	Warranty	Recom. Code
Renegade® & LP	Z-Shock technology, lightweight, sole plate, optional sandal-toe (22-28 cm)	LP: 4.7-5.6 in / 120-143 mm Tall: 6.1-7.4 in / 154-189 mm	510-515 g	365 lb / 166 kg	22-31 cm	Carbon-36 mo	L5987
Renegade® MX & LP-MX	Z-Shock technology with multi-axial rotation optional sandal-toe (22-28 cm)	LP: 5.6-6.5 in / 143-166 mm Tall: 6.7-8.1 in / 170-206 mm	850-945 g	285 lb / 130 kg	22-31 cm	Carbon-36 mo Rotator-24 mo	L5987, L5986
Thrive™ VS ²	Load activated with Dual Keel Technology, vertical shock	5.75-7.0 in / 145-175 mm	510 g	365 lb / 166 kg	22-31 cm	36 mo	L5987, L5986
Runway® & HX	User adjustable heel height foot/ankle system, push button or HX design optional sandal-toe (22-28 cm)	3.6-3.9 in / 92-98 mm	485 g	255 lb / 116 kg	22-28 cm	Carbon-36 mo Ankle-18 mo	L5981, L5990
Silhouette® VS & LP-VS	Vertical shock, multi-axial, easy to cosmetically finish	LP: 6.6-7.5 in / 168-191 mm Tall: 7.0-7.9 in / 178-201 mm	300-395 g	365 lb / 166 kg	22-31 cm	36 mo	L5987, L5986
Silhouette® & LP	Lightweight, multi-axial, easy to cosmetically finish	LP: 6.6-7.5 in / 168-191 mm Tall: 7.0-7.9 in / 178-201 mm	300-395 g	365 lb / 166 kg	22-31 cm	36 mo	L5980, L5986
Promenade®	Designed for low impact users, partial split toe	4.6-5.3 in / 116-135 mm	345 g	325 lb / 147 kg	22-31 cm	36 mo	L5981
Senator®	Sophisticated simplicity provides energy return and durability, foot shell included, Senator Categories 1-6	4.7-5.9 in / 120-149 mm	405 g	300 lb / 136 kg	22-30 cm	24 mo	L5976
Sierra®	Split toe design, longer keel, highly responsive	5.5-6.8 in / 140-172 mm	390 g	365 lb / 166 kg	22-31 cm	36 mo	L5980
Highlander®	Split toe design, favored by TF amputees	3.9-4.7 in / 98-121 mm	340 g	365 lb / 166 kg (HD up to 500 lb / 226 kg)	22-31 cm	36 mo	L5981
Pacifica® & LP	Low profile, split toe design, easy to cosmetically finish	LP: 1.9-2.3 in / 47-60 mm Tall: 2.3-2.5 in / 59-63 mm	320-360 g	365 lb / 166 kg	22-31/ 24-31 cm	36 mo	L5981
LP Symes™	Lightweight, flexible, durable with direct lamination cup	1.7-2.1 in / 43-53 mm	295 g	365 lb / 166 kg	22-31 cm	36 mo	L5981
Defender®	Flexible and durable pediatric foot for users weighing up to 100 lbs	2.4-3.0 in / 61-76 mm	107 g	100 lb / 45 kg	16-21 cm	18 mo	L5981
Nitro® Running	Custom designed for runners, tall or low profile	Low 7.25 in / 184 mm Tall 9.5 in / 241 mm	276-303 g	365 lb / 166 kg	N/A	12 mo	L5980
Slalom™ Ski	Custom made; inserts directly into ski binding	5.5 in / 140 mm	1,056 g	250 lb / 114 kg	N/A	12 mo	L5980
Freestyle™ Swim	Adjust to 70° plantar flexion; push button design	4-4.6 in / 102-117 mm	542 g	255 lb / 116 kg	22-30 cm	Carbon-36 mo Ankle-12 mo	L5981, L5990
WalkTek™	iMX™ stability, optional sandal-toe (22-28 cm)	3.4-3.8 in / 85-97 mm	630 g	300 lb / 136 kg	21-30 cm	12 mo	L5972, L5986

Silicone Socket Interfaces

Product	Key Features	Sizes	LCode	Warranty
Parydonn® Liner	3 mm silicone TT liner, cushion or locking	22, 23.5, 25, 26.5, 28, 30, 32, 34 cm	locking L5673, cushion L5679	6 mo

Dynamic Performance

Style and Flexibility

Everyday Performance

Unique by Design

Extreme Performance

K2 Level