

# Renegade A·T<sup>®</sup>

Low Profile Design Available

## All-Terrain Performance AMPlified

- The most dynamic, versatile and durable Renegade prosthetic foot
- EnduraCore™ Technology delivers 3 times the fatigue life of standard laminates
- Unique design returns up to 35% more energy than competitive alternatives
- Split-toe design delivers multi-axial function with +/- 15 degrees of inversion/eversion for enhanced stability
- Z-Shock Technology™ absorbs shear forces providing greater user comfort
- Lightweight construction – only 515 grams
- Rated to 166 kg (365 lbs), 36-month warranty

## Mobility Indicator



Renegade AT (R11)

Renegade LP-AT (R16)



# Renegade® A·T & Renegade® LP-A·T

**Product ID:** Renegade A·T (R11); Renegade LP-A·T (R16)


**Sizes:** 22-31 cm

**User Weight Rating:** 166 kg (365 lbs)

**Heel Height:** 10 mm (3/8")

**Connector:** Male pyramid

**Warranty:** Foot module 36 mo, Cosmetic foot shell 6 mo



Note: Build height measurements represent a fully assembled unit.

## Build Height & Weight

	Size (cm)	Height (in / mm)	Weight (oz / g)
Renegade A·T	22-25	6 1/8 / 154	22 / 625
	26-28	6 3/4 / 172	27 / 765
	29-31	7 1/2 / 189	33 / 945
Renegade LP-A·T	22-25	4 3/4 / 120	22 / 625
	26-28	5 1/8 / 131	27 / 760
	29-31	5 5/8 / 143	33 / 930

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

## Stiffness Category Selection

Weight (lbs)	Impact Level			Weight (kg)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9		131-147
326-365	9			148-166

## Accessories

- Pre-Cut Shin Foam: (RS1 only; two piece clamshell design)
  - 22-25 cm (AFP-00-RS101)
  - 26-28 cm (AFP-00-RS102)
  - 29-31 cm (AFP-00-RS103)
- Sole Plate: (size and category specific)
  - Universal toe: (RS1-SP)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 18"
- Stiffening Bumpers:
  - 22-25 cm: (KIT-00-116)
  - 26-28 cm: (KIT-00-117)
  - 29-31 cm: (KIT-00-118)

## Foot Module Part Number

-00-0  A  -RU  
 Product ID: R11 (Renegade A·T), R16 (Renegade LP-A·T) | Category: 1-9 | Size: 22-31

## Foot Shell Part Number

FTC-  -     4-R   
 Generation: 2-3 | Profile: M-Mid, H-High | Cap: 0-Yes, 1-No | Size: 22-31 | Skin Tone: L-Light, M-Medium, D-Dark | Side: L-Left, R-Right

See the Cosmetic Foot Shell page for more detailed ordering instructions.