


RENEGADE® LP (RS6)

The lowest profile vertical shock foot

- Less than 5" of clearance required for most commonly used sizes.
- Utilizes the same patented Z-Shock Technology™ as the standard Renegade (RS 1000)
- Less than 1/3 the weight of alternative shock foot products
- All-carbon z-shaped heel is the most durable shock design available
- Carbon fiber lay-ups are engineered to deliver a smooth gait at any speed, resulting in a comfortable walking foot that can stand up to high-activity forces
- Incorporates full-length sole plate, creating natural inversion-eversion
- Designed for K3/4 amputees
- Rated up to 365 lbs



RENEGADE® LP (RS6)

RENEGADE® LP (RS6)

Size*: 22-30 cm

*Special Order Size 31 cm (7-10 day lead time)


User Weight Rating: 365 lbs

Stiffness Category: 1-9

Heel Height: 3/8"

Connector: Integral male pyramid

Warranty: 36 months



Clearance & Weight:

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	4.4 in	337g	127g
26-28 cm	4.8 in	403g	159g
29-31 cm	5.2 in	496g	219g

Note: A Spectra sock and pyramid connector are included with purchase.

REQUIRED ACCESSORIES

Cosmetic Foot Shell:

- High-profile (Default) FTC-2H
- Mid-profile (Optional) FTC-2M
- Choice of caucasian or ebony skin tone

OPTIONAL ACCESSORIES

- Replacement Renegade Sole Plate- Part #RS6-SP
- Foot Shell Removal Tool- Part #ACC-00-103
- Additional Spectra Sock- Part #ACC-00-101xx (12", 14", 18")

Additional Stiffening Bumpers-

- KIT-00-116 40 Durometer, Sizes 22-25cm
- KIT-00-117 40 Durometer, Sizes 26-28cm
- KIT-00-118 40 Durometer, Sizes 29-31cm