


## SIERRA® PROSTHETIC FOOT (FS1)

### *Ultimate flexibility for a smoother gait*

- Angle-top design increases the length of the carbon fiber spring and resulting energy return and flexibility
- Split toe design provides stability enhancing inversion and eversion
- Full-length heel enables foot flat early in the gait cycle, providing a stable base of support and an exceptionally smooth gait
- Rated up to 365 lbs



SIERRA® (FS1)

# SIERRA® PROSTHETIC FOOT (FS1)

**Size\*:** 22-30 cm

\*Special Order Size 31 cm (7-10 day lead time)


**User Weight Rating:** 365 lbs

**Stiffness Category:** 1-9

**Connector:** Integral male pyramid

**Heel Height:** 3/8"

**Warranty:** 36 months



## Clearance & Weight:

Size	Clearance	Avg. Weights	
		Foot	Shell
22-25 cm	5.2 in	343g	127g
26-28 cm	6.0 in	404g	159g
29-31 cm	6.7 in	493g	219g

Note: A Spectra sock and pyramid connector are included with purchase.

## REQUIRED ACCESSORIES

### Cosmetic Foot Shell:

Mid-profile (Default) FTC-2M

High-profile (Optional) FTC-2H

Choice of caucasian or ebony skin tone

## OPTIONAL ACCESSORIES

- **Pre-Cut Shin Foam Two part clamshell design is size specific.**

AFP-00-FS101 For FS1, Sizes 22-25

AFP-00-FS102 For FS1, Sizes 26-28

AFP-00-FS103 For FS1, Sizes 29-31

- **Foot Shell Removal Tool- Part #ACC-00-103**

- **Additional Spectra Sock- Part #ACC-00-101 (12", 14", 18")**

- **Additional Stiffening Bumpers-**

KIT-00-110 40 Durometer, Sizes 22-25cm

KIT-00-112 40 Durometer, Sizes 26-28cm

KIT-00-114 40 Durometer, Sizes 29-31cm