

WalkTek™



Redefining K2 stability

- Offers superior stability and support with iMX Technology™
- iMX Technology delivers controlled, multi-axial movement and outstanding ground compliance to overcome small community obstacles and changes in surface angles
- A convenient, ankle-height pyramid simplifies the WalkTek fitting process by making adjustments quick and easy
- The WalkTek prosthetic foot features heavy-duty, graphite composite construction and an extremely durable, low-maintenance design, providing you peace of mind
- Rated to 136 kg (300 lbs), 12-month warranty

Mobility Indicator



WalkTek (VS3)

WalkTek™

Product ID: VS3

Sizes: 21-30cm


Sandal-toe available for sizes 22-28 cm

User Weight Rating: 136 kg (300 lbs)

Heel Height: 10 mm (3/8")

Connector: Male pyramid

Warranty: Foot module 12 mo, Cosmetic foot shell 6 mo



Note: Build height measurements represent a fully assembled unit.

Build Height & Weight

Size (cm)	Height (in / mm)	Weight (oz / g)
21-25	3 3/8 / 85	17 / 490
26-28	3 3/4 / 97	22 / 630

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Weight (lbs)	Category	Weight (kg)
100-130	1	44-59
131-180	2	60-81
181-240	3	82-109
241-300	4	100-136

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx) 12", 14", 18"

Foot Module Part Number

VS3-00-0 A -

Category: 1-4 Size: 21-30 Toe Configuration: RU-Regular toe, Universal; SR-Sandal toe, Right; SL-Sandal toe, Left

Foot Shell Part Number

FTC- - 4-

Generation: 1-3 Profile: M-Mid; H-High; S-Sandal Cap: 0-Yes; 1-No Size: 21-30 Skin Tone: L-Light; M-Medium; D-Dark Toe: R-Regular; S-Sandal Side: L-Left; R-Right

See the Cosmetic Foot Shell page for more detailed ordering instructions.