

The Promenade™

Prosthetic Foot



**Because hurdles come in all shapes
and sizes for low impact amputees**

The Promenade™ Prosthetic Foot

Designed specifically for the low impact ambulator

Freedom Innovations knows that world-class amputee athletes are the minority. We recognize that although there are no award ceremonies for everyday amputees, the hurdles they face can be as challenging as those on the track, particularly if done with the wrong prosthetic foot. Likewise, the biggest hurdle for a prosthetist is often choosing the right prosthetic foot for this patient class.

That's why Freedom Innovations developed the Promenade™ Prosthetic Foot.

The Promenade is engineered to manage the unique gait dynamics of low impact ambulators including the diabetic & vascular amputee population. Once caught in the middle between foot products that either under or over performed for this patient

category's needs, the Promenade provides an optimal balance between flexibility and stability for the utmost in comfort. With an easily deflected toe - not too hard or too soft - the product returns energy without overpowering the less mobile amputee. Users will enjoy a more natural and efficient stride, allowing them to increase their daily activity.

While taking a stroll on the beach may not be an Olympic event, in the real world everyday hurdles can be just as challenging. Fortunately there is one prosthetic foot that is engineered with you and your patient in mind - the Promenade™ Prosthetic Foot.



Near effortless roll-over optimized for low impact amputees

A dynamic heel and easily deflected toe provide near effortless roll-over, helping to reduce gait deviations often presented by low impact ambulators. As a result, amputees save energy so they can do more and feel better at the end of the day.



Partial split toe

The split toe design offers controlled inversion/eversion for stability on uneven surfaces. Users will gain the confidence needed for a variety of walking conditions.



Full Length Sole Plate

The full length carbon sole plate provides increased energy storage and return and a stable base of support, promoting confidence among users. A high strength wrap attaches the sole plate to the foot module, eliminating the stress points and 'dead spots' associated with bolted, heel only products. Users experience an uninterrupted, continuously smooth roll-over and increased comfort throughout their daily activities.

High User Weight Rating

The robust, 325 lb. user weight rating makes the Promenade appropriate for a broader range of users compared to competitive alternatives.



SIMPLE

Simplified Fitting Process

An integral male pyramid and just 5 weight categories from which to choose make for an easy fitting process, saving you time and money so you can move on to what's next.

Warranty, Trial Period

A 36 month warranty and 60 day trial period protect your investment and provide peace of mind over the duration of the product's life.



CATEGORY SELECTION CHART

Weight (lbs)	Impact Level	Weight (kg)
	Low	
100-126	1	45-57
127-160	2	58-72
161-203	3	73-92
204-258	4	93-117
259-325	5	118-147

PRODUCT CLEARANCE & WEIGHT

Size (cm)	Clearance	Weight (kg)	
		Foot	Shell
22-25	4.8 in	280 g	127 g
26-28	5.2 in	330 g	159 g
29-31	5.6 in	375 g	219 g

**Special order size 31cm*



SPECIFICATIONS

Size offering - 22-30 cm

**Special order size 31 cm*

Average clearance - 5.2 in.

Average weight - 330 g

Carbon fiber foot warranty - 36 months

Foot shell warranty - 6 months

Suggested LCode - L5981



30 Fairbanks • Ste. 114 • Irvine, CA 92618 • toll free: 888-818-6777 • ph: 949-672-0032
www.freedom-innovations.com

Also distributed by SPS, call (800) 767-7776